



Exceeding Expectation

Triangle Surgical is the UK's leading supplier of Liposuction equipment and instrumentation, we are also very proud of the trading partnerships forged in our ever expanding export market.

A family company that believes strongly in the quality of our products and the longevity of our relationships, our aim is to exceed your expectations and you can trust our advice as we have over 100 years of combined experience in this speciality alone.



The outstanding quality of the CETTA disposable product range is immediately apparent; high grade stainless steel, reassuring weight, excellent balance and superb finish. We strongly believe CETTA is not only the best disposable product available worldwide but is also far superior to the majority of re-usable brands. For Cetta products please refer to the product codes in blue.

The reusable Primo range is as versatile as you need it to be, we will create any bespoke size and style upon request. The range shown in this catalogue is a guide to the most popular products we currently provide, many are regular stock items for fast delivery. Durable, high quality Stainless Steel shafts, reassuringly weighty and balanced Aluminium handles. For Primo products please refer to the product codes in gold.

FLUID INFUSION

Infusion needles

High grade stainless steel. 6 spaced holes in a spiral around the tip. For fluid infusion using leur lock syringes, handles or giving set tubing.



CETTA RANGE	PRIMO RANGE		
Product Code	Product Code	Diameter (mm)	Length (cm)
CETTA INF/N 215	PRIMO INF/N 215	2	15
CETTA INF/N 223	PRIMO INF/N 223	2	23
CETTA INF/N 230	PRIMO INF/N 230	2	30
CETTA INF/N 315	PRIMO INF/N 315	3	15
CETTA INF/N 323	PRIMO INF/N 323	3	23
CETTA INF/N 330	PRIMO INF/N 330	3	30

Product Code	Product Code	Diameter (mm)	Length (cm)
CETTA INF/N1.623CRV	PRIMO INF/N1.623CRV	1.6	23

Infusion handle non valve

Comfortable handle to link tubing to needle, lack of valve reduces leakage. Single shaft allows infiltration and suction.

Product Code

CETTA INF/HAN W/O

PRIMO INF/HAN W/O



Infusion handle with Valve

This infusion handle, with on/off button, offers maximum control for infusion during the tumescent technique.

Product Code

PRIMO HUN/HAN

02 |



POWERED SUCTION

One piece cannulae

High grade Stainless Steel shaft, for aspiration with powered suction.



Primo handle: Aluminium



Cetta handle: POM

CETTA RANGE	PRIMO RANGE		
Product Code	Product Code	Diameter (mm)	Length (cm
CETTA 210MM	PRIMO 210MM	2	10
CETTA 315MM	PRIMO 315MM	3	15
CETTA 325MM	PRIMO 325MM	3	25
CETTA 430MM	PRIMO 430MM	4	30
CETTA 530MM	PRIMO 530MM	5	30

Product Code	Product Code	Diameter (mm)	Length (cm)
CETTA 325/30holeM	PRIMO 325/30holeM	3	25

Product Code	Product Code	Diameter (mm)	Length (cm)
CETTA 325LM	PRIMO 325LM	3	25
CETTA 430LM	PRIMO 430LM	4	30
CETTA 530LM	PRIMO 530LM	5	30

0
0
0
)

r (mm) Length (cm
10
15

Product Code	Product Code	Diameter (mm)	Length (cm
CETTA 215OHM	PRIMO 2150HM	2	15
CETTA 3200HM	PRIMO 3200HM	3	20
CETTA 4300HM	PRIMO 4300HM	4	30

Mercedes

Three holes spaced evenly around the tip. The most popular tip style.



30 Hole harvest cannula

Harvest fat in smaller parcels, reduces blocking.



Lantern

Chinese lantern shaped tip, useful for fibrous areas and popular in Gynecomastia procedures.



The Scream

Two holes above a single hole, a cannula lookalike of the famous Munch painting.



Submental

One hole, flat backed for areas closer to the skin.



One hole

A single downward facing aperture.

cing aperture.



TOOMEY SYRINGE SUCTION



Length (cm)

Toomey Cannulae

CETTA RANGE

Product Code

CETTA 210MTS

High grade Stainless steel shaft,

for aspiration with 60cc Toomey syringes.

PRIMO RANGE

Product Code

PRIMO 210MTS

MICRO HARVEST

Tissue harvest using luer lock syringes. All Triangle Surgical luer lock bases have a thumb platform for comfort and orientation.

Open harvest

Guarded fully open tip, designed specifically for tissue harvest.

CETTA RANGE	PRIMO RANGE		
Product Code	Product Code	Diameter (mm)	Length (cm)
CETTA MIC315/HARV	PRIMO MIC315/HAR	V 3	15
CETTA MIC322/HARV	PRIMO MIC322/HARY	V 3	22

Product Code



Primo handle: Aluminium

Mercedes tip

Three holes spaced evenly around the tip. The most popular tip style.



30 Hole harvest cannula

Harvest fat in smaller parcels, reduces blocking.




Submental

One hole, flat backed for areas closer to the skin.



Open harvest

Guarded fully open tip, designed specifically for tissue harvest.



Diameter (mm)

Product Code	Product Code	Diameter (mm)	Length (cm)
CETTA 2.515/30holeTS	PRIMO2.515/30holeTS	2.5	15
CETTA 320/30holeTS	PRIMO320/30holeTS	3	20

Product Code	Product Code	Diameter (mm)	Length (cm)
CETTA 210SMTS	PRIMO 210SMTS	2	10
CETTA 315SMTS	PRIMO 315SMTS	3	15

Product Code	Product Code	Diameter (mm)	Length (cm)
CETTA 318HARV/TS	PRIMO 318HARV/TS	3	18

Guarded fully open tip, c



All 60cc locks also fit 50cc BD

60cc syringe lock - single arm

Product Code	Product Code	Syringe Size(co
CETTA SNAP10	PRIMO SNAP10	10cc
CETTA SNAP20	PRIMO SNAP20	20cc
CETTA SNAP60	PRIMO SNAP60	60cc

60cc syringe lock - double arm

Fixes to syringe base for added security.

Product Code	Product Code	Syringe Size(cc)
CETTA DOUBSNAP60	PRIMO DOUBSNAP60	60cc

04

60cc syringe lock - 2 step

Product Code	Product Code	Syringe Size(co
CETTA2STEP	PRIMO2STEP	60cc

Square syringe locks Product Code Syringe Size(cc)

roduct Code	Syringe Size(cc)
RIMOSQUARE10	10cc
RIMOSQUARE30	30cc
RIMOSQUARE60	60cc



* *



600000

30 Hole harvest cannula

Harvest fat in micro droplets for fine graft re-injection

at in micro dronlets for tine graft re-injection	PRIMO 2.515LL/30hole PRIMO 315LL/30hole	2.5 3	15 15

Product Code

Submental

One hole, flat backed for areas closer to the skin.

Product Code	Product Code	Diameter (mm)	Length (cm)
CETTA 322S/LL	PRIMO 322S/LL	3	22

Mercedes tip

Three holes spaced evenly around the tip. The most popular tip style.

Product Code	Product Code	Diameter (mm)	Length (cm)
CETTA 215MLL	PRIMO 215MLL	2	15
CETTA 315MLL	PRIMO 315MLL	3	15
CETTA 322MLI	PRIMO 322MLI	3	22



CELL SELECT CANNULAE RANGE

Level 1 - Gentle Rasp

		-			
Product Code	Product Code	Holes	Diameter (mm)	Length (cm)	
CETTA 2515G	PRIMO 2515G	21	2.5	15	
CETTA 320G	PRIMO 320G	30	3	20	

Level 2 - Medium Rasp

Product Code	Product Code	Holes	Diameter (mm)	Length (cm)
CETTA 2515M	PRIMO 2515M	21	2.5	15
CETTA 320M	PRIMO 320M	30	3	20

Level 3 - Aggressive Rasp

			1-		
Product Code	Product Code	Holes	Diameter (mm)	Length (cm)	
CETTA 2515A	PRIMO 2515A	21	2.5	15	
CETTA 320A	PRIMO 320A	30	3	20	

The ingenious Cell Select micro harvest cannula range is specifically designed to capture smaller aliquots of healthy fat cells via numerous 1mm micro-ports. The harvest of smaller cell packages results in easier, smoother re-injection through finer injection cannula. The spacing and increased surface area of the apertures per cannula tip results in an effective and constant collection speed, and greatly reduces blocking.











MICRO INFILTRATION CANNULAE

High grade Stainless Steel for re-injection using Luer lock syringes.

CETTA RANGE	PRIMO	RANGE
-------------	-------	-------

Product Code	Product Code	Diameter (mm)	Length (cn
CETTA MIC16/A7	PRIMO MIC16/A7	16G	7
CETTA MIC16/A9	PRIMO MIC16/A9	16G	9
CETTA MIC16/B7	PRIMO MIC16/B7	16G	7
CETTA MIC16/B9	PRIMO MIC16/B9	16G	9
CETTA MIC16/B12	PRIMO MIC 16/B12	16G	12
CETTA MIC16/C7	PRIMO MIC16/C7	16G	7
CETTA MIC16/C9	PRIMO MIC16/C9	16G	9
CETTA MIC18/A9	PRIMO MIC18/A9	18G	9
CETTA MIC18/B9	PRIMO MIC18/B9	18G	9



Curved Cannula

Product Code	Product Code	Diameter (mm)	Length (cm
CETTA CURV/18B7	PRIMO CURV/18B7	18G	7
CETTA CURV/16B12	PRIMO CURV/16B12	16G	12
CETTA CURV/16B15	PRIMO CURV/16B15	16G	15

Ultra Fine

Ultra fine facial needle with unblocking Stylet

Product Code	Product Code	Diameter (mm)	Length (cm
CETTA 235FN/LL	PRIMO 235FN/LL	23G (0.7mm)	5
CETTA 217FN/LL	PRIMO 217FN/LL	21G (0.8mm)	7
CETTA 2185FN/LL	PRIMO 2185FN/LL	21G (0.8mm)	8.5
CETTA 15FN/LL	PRIMO 15FN/LL	20G (1.0mm)	5
CETTA 207FN/LL	PRIMO 207FN/LL	20G (1.0mm)	7
CETTA 189FN/LL	PRIMO 189FN/LL	18G (1.2mm)	9



B tip

A tip


C tip

Micro Drop

Curved obtuse tip for accurate placement of micro droplets.

Product Code	Product Code	Diameter (mm)	Length (cm)	
CETTA 1285MD	PRIMO 1285MD	1.2	8.5	





DISSECTORS

V-Dissector

Undermining fibrous tissue, widens at the tip and ideal for breast surgery.

CETTA RANGE	PRIMO RANGE		
Product Code	Product Code	Diameter (mm)	Length (cm
CETTA VDISS167	PRIMO VDISS167	16G	7
CETTA VDISS169	PRIMO VDISS169	16G	9

W-Dissector

For undermining small areas of fibrous tissue, less aggressive and finer than V Dissector.

Product Code	Product Code	Diameter (mm)	Length (cm)
CETTA WDISS	PRIMO WDISS	16G	8.5



HAY Tunnel Dissector

For blunt dissection tunnelling when a V-dissector may be considered aggressive i.e. around Implants. Slight upward curve to tip in order to inhibit direct forceful contact, comfortable handle with orientation marking.

Product Code	Product Code	Diameter (mm)	Length (cm)
CETTA HTD	PRIMO HTD	2	12

ADAPTORS

Toomey to luer lock adaptor

Transfer aspirate from 60cc Toomey to any size luer lock syringe.

Produc	t Code	Product Code						
CETTA	MIC/LL/S	PRIMO MIC/LL/S						
	Skila	lanhadaripadaripada	II.		2	1 1	02	22

CETTA MIC/LL/LL PRIMO MIC/LL/LL

Luer lock to luer lock adaptor Transfer between any size luer lock syringes Product Code

Toomey syringe to luer lock needle adaptor

Product Code	Product Code
CETTA MIC/LL/N	PRIMO MIC/LL/



ACCESSORIES

Trigger injection gun 10cc/5cc

Trigger mechanism controlled injection from 10cc BD syringe or 5cc BD syringe (adaptor included), eases the physical pressures of re-injection and controls accidental over filling.

Product Code		Capacity (cc)		
	PRIMO TG105	5. 10		



Syringe rack


Autoclavable steel, holds 60 10cc syringes or 12 60cc syringes

Product Code PRIMO MIC/10CCR PRIMO MIC/60CCR

The golden gun 10cc

Enjoy being Scaramanga with your very own golden gun. Thumb controller injection from 10cc BD syringe, eases the physical pressures of re-injection and controls accidental over filling.

Product Code	Capacity (co
PRIMO GG10	10cc









Packing and Certification:

Double packed, ETO sterilised 5 year shelf life, four traceability stickers included.





Master Indemnity Registration number is IFA814

Purchasing Systems Include





Quality Validation





All of our products meet the required CE certification. CETTA is grade 2 sterile product certified.

Tel: 0044 (0)1494 876 990 Fax: 0044 (0)1494 876 991 Email: info@trianglesurgical.uk

6-10 Church Farm Courtyard High Street Chalfont St Giles Buckinghamshire HP8 4QH UK

www.trianglesurgical.uk

