

REJUVENATION SOLUTIONS

Micro Fat • Nano Fat (SVT)

State of the Art

The exclusive or combined use of micro fat and nano fat (SVT) in regenerative and rejuvenation therapies is taking our industry by storm.

The plastic surgeon's access to autologous tissue harvest has resulted in the speciality leading the way in the development of the skills and equipment necessary to enact cellular rejuvenation.

Here we have assembled the very latest state of the art products which combined or used in their own right provide the ultimate tools for use within this field of medicine.

Rejuvenate • Regenerate • Repair

- Face
- Hands
- Hair
- Neck
- Décollotage
- Genital

Raising the Standard




TRIANGLE
SURGICAL LTD

State of the art Rejuvenation solutions - Raising the Standard

CANNULA

Thrive - Rejuvenation cannula. Beautiful instruments designed to offer the ideal harvest and re-injection requirements of rejuvenation surgery.

- Harvest Cannula - Standardising and sizing the aspirate via 20 gently bevelled 1mm micro ports.
- Re-injection Cannula - Ultra fine non-flex stainless steel, blunt tipped. 3 Sizes to accurately, comfortably and safely deliver Macro, Micro and micronized Nano fat cells.
- Luer to Luer connector.
- Lockable harvest syringe



Thrive
Specialist Cannula



- ✓ Single use
- ✓ Individually packed
- ✓ Superior quality
- ✓ Non-flex safety
- ✓ No plastic hubs
- ✓ Ergonomic design
- ✓ World-wide popularity since 2011

BEAUTY-STEM

BEAUTY-STEM is an ingenious closed-circuit, single use kit for adipose tissue purification without the need for any centrifuge and with minimal manipulation. With its simple Lipo In and Lipo Out valve system and continuous wash process the BEAUTY-STEM is really simple to use and requires only one operator with minimal training.

The unique filtration system ensures the purification process eliminates all unwanted contaminants and preserves the biological properties of the tissue, reducing cellular stress and maximising its regenerative potential.



BEAUTY-STEM™
SUPERIOR ADIPOSE TISSUE



- ✓ Ready to go, single use kit
- ✓ Sterile, closed circuit environment
- ✓ Superior adipose tissue
- ✓ Simple to use, minimal training required
- ✓ Continual wash enables processing of up to 400ml of adipose tissue in 10 minutes
- ✓ Easily retrievable graft



**TIME
EFFICIENT**



**MINIMAL
MANIPULATION**




EASY



State of the art Rejuvenation solutions - Raising the Standard

FAT LOCK ONE WAY VALVE

Fully closed tubing and adaptor system for collection of graft with luer lock syringes. No need, even, to remove the cannula from the harvest site to re-load your suction.



- ✓ Unique valve & tubing system kit
- ✓ Closed system syringe harvest
- ✓ Saves valuable time
- ✓ Use with lockable harvest syringe (available separately)



Disposable self-locking syringe. Twist locks at 5cc increments along the barrel.

EMULSFAT SVT PREPARATION

The Emulsfat is an innovative and ergonomically designed device for the emulsification of adipose tissue and then the filtration of this emulsified tissue to produce a high yield of high quality SVT, ready to be re-injected.

The ergonomic 'Y' shape adopts a natural position for the hands providing effortless handling for the practitioner and promoting an easy and consistent technique.

Emulsfat is a closed system, requiring only between 10 - 20 passes for the emulsification of the adipose tissue, depending on the required SVT fluidity.



- ✓ Simple, easy, fast SVT preparation
- ✓ Closed system process
- ✓ Unique ergonomic design
- ✓ 80% cell viability
- ✓ Over 90% SVT yield
- ✓ Only 10-20 passes required

ADDITIONS

Choose from our range of top quality consumable items specially selected to provide your one stop shop service for a complete rejuvenation package offering. All items are available for speedy delivery and at a great price.

Stab incision needle, 21G x 50mm	✓
Meso needle 27g 4mm	✓
Meso multi injector 27g 4mm x 3	✓
1ml BD syringe	✓
10ml BD syringe	✓
20ml BD syringe	✓



ORDER CODES	DESCRIPTION
CANNULA	
THRIVE-001	Thrive specialist rejuvenation cannula set 1 x CETTA-2515/20HOLE/LL - harvest cannula 1 x CETTA-235FN/LL - ultra fine injection cannula with unblocking stylet 1 x CETTA-217FN/LL - ultra fine injection cannula with unblocking stylet 1 x CETTA-189FN/LL - ultra fine injection cannula with unblocking stylet 1 x CETTA-MIC/LL/LL - luer to luer connector 1 x SYRINGE-30ML-LOCK - lockable harvest syringe
BEAUTY-STEM	
BETKIT-004	BEAUTY-STEM adipose tissue purification system kit
FAT LOCK SYSTEM	
PURE-FLS	Fat Lock system-single use one-way valve and tubing set
SYRINGE-30ML-LOCK	Lockable harvest syringe. Twist locks at 5cc increments along the barrel.
EMULSFAT SVT PREP	
CGHAPI-Y	Emulsfat adipose tissue emulsifier and filtration system for SVT preparation
ADDITIONS	
SIN-2150	Stab incision needle, 21G 50mm
AAL4	Meso needle, 27g, 4mm
MRLA273	Meso multi injector, 27g, 4mm
SYRINGE-1ML-BD	1ml BD luer lock tip syringe, box of 100
SYRINGE-10ML-BD	10ml BD luer lock tip syringe
SYRINGE-20ML-BD	20ml BD luer lock tip syringe

Other speciality cooperations include:

Sports Medicine
Dermatology
Orthopedics
Gastroenterology
Neurology
Urogenital

Triangle Surgical Ltd
6-10 Church Farm Courtyard
Chalfont St Giles
Buckinghamshire
HP8 4QH, UK

Tel: +44 (0) 1494 876990
e: info@trianglesurgical.uk



TRIANGLE
SURGICAL LTD